

1


2

3

4

## ISEF121 probe

## Introduction

The probe is a combination fluoride selective electrode with a non-refillable reference and built-in temperature sensor.

**Warning**

**⚠ Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.**

**User manual**

A user manual with all instructions is on the manufacturer's website.

**Collect the following items:**

Compatible meter, Fluoride standard solutions and ionic strength adjusters; beakers, 50 mL; stirrer plate; DI water; stir bars; lint-free cloth.

**1 Connect the probe****2 Calibration**

*Air bubbles under the probe tip when submerged can cause slow stabilization or error in measurement. If bubbles are present, gently shake the probe until bubbles are removed.*

**3 Measurement (direct method)**

Temperatures of calibration standards and samples should be kept within  $\pm 2$  °C of each other for optimal results.

**4 Storage****Technical information**

Range	0.01 mg/L ( $1 \times 10^{-6}$ ) to 19,000 mg/L (1 M) F <sup>-</sup>
pH Range	pH 4 to 8, adjusted to 5.0 to 5.5 by ISA
Operating temperature	5 to 50 °C (41 to 122 °F)
Minimum sample volume	25 mL
Ionic strength adjusters	Hach fluoride ISA (1 powder pillow per 25 mL standard or sample)
Default standard solutions	0.5, 1 and 2 mg/L Hach fluoride standard solutions

**Maintenance and storage**

The probe can be stored dry. To protect the sensing element, rinse with DI water and blot dry with a lint-free cloth. Install the sensor protection cap.

**Warranty**

1 year for manufacturing faults. This warranty does not cover improper use or wear.

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