



# CiDehol<sup>®</sup> 70

## 70% IPA Solution



**CiDehol<sup>®</sup> 70** is a 70% (v/v) Isopropyl alcohol solution filtered to 0.2 microns. **CiDehol<sup>®</sup> 70** can be used to clean surfaces in microbiology labs and tissue culture labs. It is also effective in cleaning surfaces and process equipment in pharmaceutical, biotechnology, and medical device facilities. IPA solution evaporates away completely so no rinse is required.

**Use Protocol:** Use in well-ventilated area. Spray or wipe solution on surface to be cleaned until it is completely covered. Allow to air dry. When using trigger spray, point away from eyes and skin. Hold bottle six to eight inches from surface to be treated.

### *Precautionary Statements*

Flammable - do not come in contact with flame or extreme heat. Avoid eye and skin contact and ingestion.

### *INGREDIENTS*

70 % Isopropyl Alcohol, USP

30% Purified Water, USP

**SHELF LIFE**  
**2 YEARS**

Size	Case Qty	Decon Cat #
16 oz.	12	8416
16 oz. Color Label	12	8417
32oz.	12	8432
1 Gallon	4	8401
5 Gallon	-	8406
55 Gallon	-	8405

Rev. 07/2021